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## 11 Anhang Bibliografie

1. Ardan NI Jr, Janes JM, Herrick JF. Ultrasonic energy and defects in bone. *J Bone Joint Surg Am* 1957; 39:394-402
2. Ay S, Dogan SK, Evcik D, Baser OC. Comparison the efficacy of phonophoresis and ultrasound therapy in myofascial pain syndrome. *Rheumatol Int.* 2011; 31(9):1203-08
3. Bakhtiary AH, Rashidy-Pour A. Ultrasound and laser therapy in the treatment of carpal tunnel syndrome. *Aust J Physiother.* 2004; 50(3):147-51
4. Bashardoust Tajali SB, Houghton P, MacDermid JC, Grewal R. Effects of low-intensity pulsed ultrasound therapy on fracture healing: a systematic review and meta-analysis. *Am J Phys Med Rehabil* 2012; 91(4):349-67
5. Basso O, Pike JM. The effect of low frequency, long-wave ultrasound therapy on joint mobility and rehabilitation after wrist fracture. *J Hand Surg Br* 1998; 23:136-9
6. Baysal O, Altay Z, Ozcan C, Ertem K, Yologlu S, Kayhan A. Comparison of three conservative treatment protocols in carpal tunnel syndrome. *Int J Clin Pract.* 2006; 60(7):820-8
7. Bender T, Gidófalvi E. the effect of ultrasonic therapy in rheumatoid arthritis of the temporomandibular joint. *Fogorv Sz.* 1991; 84(8):229-32
8. Bilgici A, Ulusoy H, Kuru O, Canturk F. The comparison of ultrasound treatment and local steroid injection plus splinting in the carpal tunnel syndrome: a randomized controlled trial. *Bratisl Lek Listy.* 2010; 111(12):659-65
9. Brand JC, et al. Does pulsed low intensity ultrasound allow early return to normal activities when treating stress fractures? *Iowa Orthopaedic Journal* 1999; 19:26-30
10. Busse JW, et al. The effect of low-intensity pulsed ultrasound therapy on time to fracture healing: A meta-analysis. *Canadian Medical Association Journal* 2002; 166:437-441
11. Cambier D, D'Herde K, Witvrouw E, Beck M, Soenens S, Vanderstraeten G. Therapeutic ultrasound temperature increase at different depths by different modes in a human cadaver. *J Rehabil Med.* 2001; 33(5):212-5
12. Chang CJ, Hsu SH, Lin FT, Chang H, Chang CS. Low-intensity-ultrasound-accelerated nerve regeneration using cell-seeded poly(D,L-lactic acid-co-glycolic acid) conduits: an in vivo and in vitro study. *J Biomed Mater Res B Appl Biomater.* 2005; 75(1):99-107
13. Cook SD, Ryaby JP, McCabe J, Frey JJ, Heckman JD, Kristiansen TK. Acceleration of tibia and distal radius fracture healing in patients who smoke. *Clin Orthop* 1997; 337:198-207
14. Crisci AR, Ferreira AL. Low-intensity pulsed ultrasound accelerates the regeneration of the sciatic nerve after neurotomy in rats. *Ultrasound Med Biol.* 2002; 28(10):1335-41
15. Dakowicz A, Latosiewicz R. The value of iontophoresis combined with ultrasound in patients with the carpal tunnel syndrome. *Rocz Akad Med Bialymst.* 2005; 50 (Suppl 1):196-8
16. Demmink JH, Helders PJM, Hobæk H, Enwemeka C. The variation of heating depth with therapeutic ultrasound frequency in physiotherapy. *Ultrasound Med Biol.* 2003; 29(1):113-118
17. Denz S, Topuz O, Atalay NS, Sarsan A et al. Comparison of the effectiveness of pulsed and continuous diclophenac phonophoresis in treatment of knee osteoarthritis. *J Phys Ther Sci.* 2009;

18. Dincer U, Cakar E, Kiralp MZ, Kilac H, Dursun H. The effectiveness of conservative treatments of carpal tunnel syndrome: splinting, ultrasound, and low-level laser therapies. *Photomed Laser Surg.* 2009; 27(1):119-25
19. Duarte LR. The stimulation of bone growth by ultrasound. *Arch Orthop Trauma Surg* 1983;101: 153-9
20. Dudda M, Hauser J, Muhr G, Esenwein SA. Low-intensity pulsed ultrasound as a useful adjuvant during distraction osteogenesis: a prospective, randomized controlled trial. *J Trauma.* 2011 71(5):1376-80
21. Dündar Ü, Solak Ö, Şamlı F, Kavuncu V. Effectiveness of Ultrasound Therapy in Cervical Myofascial Pain Syndrome: A Double Blind, Placebo-Controlled Study. *Turk J Rheumatol.* 2010; 25(3): 110-5
22. Dyson M, Brookes M. Stimulation of bone repair by ultrasound. *Ultrasound Med Biol* 1983;Suppl 2:61-6
23. Dyson M. Mechanisms involved in therapeutic ultrasound. *Physiotherapy.* 1987; 73:116-20
24. Dyson M. Therapeutic applications of ultrasound. In: Nyborg WL Ziskin MC. *Biological effects of ultrasound.* New York: Churchill Livingstone 1985: 121-133
25. Ebenbichler GR, Erdogmus CB, Resch KL, Funovics MA et al. Ultrasound therapy for calcific tendinitis of the shoulder. *The New England Journal of Medicine* 1999; 340(20):1533-8
26. Ebenbichler GR, Resch KL, Nicolakis P, Wiesinger GF, Uhl F, Ghanem AH, Fialka V. Ultrasound treatment for treating the carpal tunnel syndrome: randomized "sham" controlled trial. *BMJ.* 1998; 316(7133):731-5
27. Ebrahimi S, Abbasnia K, Motealleh A, Kooroshfard N, Kamali F, Ghaffarinezhad F. Effect of lidocaine phonophoresis on sensory blockade: pulsed or continuous mode of therapeutic ultrasound? *Physiotherapy.* 2012 Mar; 98(1):57-63
28. El-Mowafi H, Mohsen M. The effect of low-intensity pulsed ultrasound on callus maturation in tibial distraction osteogenesis. *Int Orthop.* 2005;29(2):121-4
29. Emami A, Petré-Mallmin M, Larsson S. No effect of low-intensity ultrasound on healing time of intramedullary fixed tibial fractures. *J Orthop Trauma* 1999;13:252-7
30. Ennis WJ, Valdes W, Gainer M, Meneses P. Evaluation of clinical effectiveness of MIST ultrasound therapy for the healing of chronic wounds. *Adv Skin Wound Care.* 2006 Oct;19(8):437-46
31. Enwemeka CS. The effects of therapeutic ultrasound on tendon healing. A biomechanical study. *Am J Phys Med Rehabil.* 1989 Dec;68(6):283-7
32. Erickson RI. Ultrasound: a usefull adjunct in temporomandibular therapy. *Oral Surg Oral Med Oral Pathol.* 1964;18:176-9
33. Everbach EC, Francis CW. Cavitation mechanisms in ultrasound-accelerated thrombolysis at 1 MHz. *Ultrasound Med Biol.* 2000 Sep;26(7):1153-60
34. Frizzell, L. A. and F. Dunn (1982). *Biophysics of ultrasound. Therapeutic Heat and Cold.* J. Lehmann. Baltimore, Williams & Wilkins

35. Fu SC, Shum WT, Heung LK, MD, Margaret Wan-Nar Wong, MBBS, Ling Qin, PhD, and Kai-Ming Chan, MD. Low-Intensity Pulsed Ultrasound on Tendon Healing A Study of the Effect of Treatment Duration and Treatment Initiation. *Am J Sports Med.* 2008; 636(9):1742-9
36. Giombini A, Di Cesare A, Casciello G, Sorrenti D, Dragoni S, Gabriele P. Hyperthermia at 434 MHz in the treatment of overuse sport tendinopathies: a randomised controlled clinical trial. *Int J Sports Med.* 2002 Apr;23(3):207-11
37. Gray RJ, Quayle AA, Hall CA, Schofield MA. Physiotherapy in the treatment of temporomandibular joint disorders: a comparative study of four treatment methods. *Br Dent J.* 1994;176(7):257-61
38. Grubisic F, Grazio S, Jajic T, Nemcic T. Therapeutic ultrasound in chronic low back pain treatment. *Journal Article: Reumatizam* 2006; 53(1):18-21
39. Hadjiargyrou M, McLeod K, Ryaby JP, Rubin C. Enhancement of fracture healing by low intensity ultrasound. *Clin Orthop* 1998; 355(Suppl):s216-29
40. Hanneman PF et al. The clinical and radiological outcome of pulsed electromagnetic field treatment for acute scaphoid fractures: a double-blind placebo-controlled multicentre trial. *J. Bone Joint Surg Br.* 2012;94:1403-8
41. Heckman JD, Ryaby JP, McCabe J, Frey JJ, Kilcoyne RF. Acceleration of tibial fracture-healing by non-invasive, low-intensity pulsed ultrasound. *Journal of Bone and Joint Surgery* 1994; 76A:26-34
42. Heckman JD, Sarasohn-Kahn J. The economics of treating tibia fractures. The cost of delayed unions. *Bull Hosp Joint Dis* 1997;56:63-72
43. Hough DR, Schneider HE, Rectus sheath hematoma with pain relieve after ultrasound: case report. *Mil Med* 1983;148(11):885-6
44. Huang MH, Lin YS, Lee CL, Yang RC. Use of ultrasound to increase effectiveness of isokinetic exercise for knee osteoarthritis. *Arch Phys Med Rehabil.* 2005;86(8):1545-51
45. Huisstede BM, Hoogvliet P, Randsdorp MS, Glerum S, van Middelkoop M, Koes BW. Carpal tunnel syndrome. Part I: effectiveness of nonsurgical treatments--a systematic review. *Arch Phys Med Rehabil.* 2010;91(7):981-1004
46. Hurst LC, Badalamente MA. Non operative treatment of Dupuytren's disease. *Hand Clin.* 1999;15(1):97-107
47. Jensen JE. Stress fracture in the world class athlete: a case study. *Med Sci Sports Excerc* 1998; 30: 783-787
48. Jeremias Júnior SL, Camanho GL, Bassit AC, Forgas A, Ingham SJ, Abdalla RJ. Low-intensity pulsed ultrasound accelerates healing in rat calcaneus tendon injuries. *J Orthop Sports Phys Ther.* 2011;41(7):526-31
49. Josza L, Kannus P. Human tendons. *Anatomy, physiology and pathology.* Champaign, IL, 1997
50. Kahn J, editor. *Ultrasound.* In: *Principles and practice of electrotherapy.* 4th ed. New York: Churchill Livingstone; 2000. p. 49-68
51. Khanna A, Nelmes RTC, Gougoulas N, Maffulli N et al. the effects of LIPUS on soft tissue healing: a review of the literature. *British Medical Bulletin* 2009;89:169-82

52. Kim TY, Jung DI, Kim YI, Yang JH, Shin SC. Anesthetic effects of lidocaine hydrochloride gel using low frequency ultrasound of 0.5 MHz. *J Pharm Pharm Sci.* 2007; 10(1):1-8
53. Klaiman MD, Shrader JA, Danoff JV et al. Phonophoresis versus ultrasound in the treatment of common musculoskeletal conditions. *Med Sci Sports Exerc.* 1998, 30:9, 1349-55
54. Korstjens CM, Nolte PA, Klein-Nulend J, Albers GHR, Burger EH. De invloed van ultrageluid van lage intensiteit op het bot. Een overzicht. *Ned Tijdschr Osteoporose en andere Botziekten* 2002; 6(2): 27-31
55. Korstjens CM, Nolte PA, Klein-Nulend J, Albers GHR, Burger EH. Effecten van lage-intensiteit ultrageluid op bot. *Perspectieven voor de tandheelkunde?* *Ned Tijdschr Tandheelkd* 2002; 109: 485-9
56. Kozanoglu E, Basaran S, Guzel R, Guler-Uysal F. Short term efficacy of ibuprofen phonophoresis versus continuous ultrasound therapy in knee osteoarthritis. *Swiss Med Wkly* 2003;133:333–8
57. Kristiansen TK, Ryaby JP, McCabe J, Frey JJ, Roe LR. Accelerated healing of distal radius fractures with the use of specific, low-intensity ultrasound. *Journal of Bone and Joint Surgery* 1997; 79A: 961-73
58. Kwolek A, Zwolinska J. Immediate and long-term effects of selected physiotherapy methods in patients with carpal tunnel syndrome. *Ortop Traumatol Rehabil.* 2011; 13(6):555-64 (Article in English, Polish)
59. Leung MC, Ng GY, Yip KK. Effect of ultrasound on acute inflammation of transected medial collateral ligaments. *Arch Phys Med Rehabil.* 2004;85(6):963-6
60. Loyola-Sánchez A, Richardson J, Beattie KA, Otero-Fuentes C, Adachi JD, MacIntyre NJ. Effect of low-intensity pulsed ultrasound on the cartilage repair in people with mild to moderate knee osteoarthritis: a double-blinded, randomized, placebo-controlled pilot study. *Arch Phys Med Rehabil.* 2012; 93(1):35-42
61. Loyola-Sánchez A, Richardson J, MacIntyre NJ. Efficacy of ultrasound therapy for the management of knee osteoarthritis: a systematic review with meta-analysis. *Osteoarthr. Cartil.* 2010;18(9): 1117-26
62. Majlesi J, Unalan H. Effect of treatment on trigger points. *Curr Pain Headache Rep.* 2010;14(5): 353-60
63. Majlesi J, Unalan H. High-power pain threshold ultrasound technique in the treatment of active myofascial trigger points: a randomized, double-blind, case-control study. *Arch Phys Med Rehabil.* 2004; 85(5):833-6
64. Malizos KN, et al., Low-intensity pulsed ultrasound for bone healing: an overview. *Injury.* 2006;37 (Suppl 1):S56-62
65. Markert, CD, Merrick MA, et al. Nonthermal ultrasound and exercise in skeletal muscle regeneration. *Arch Phys Med Rehabil.* 2005; 86(7):1304-10
66. Mayr E, Frankel V, Rüter A. Ultrasound: An alternative healing method for nonunions? *Archives of Orthopaedic and Trauma Surgery* 2000; 120:1-8
67. Mayr E, Rudzki MM, Rudzki M, Borchardt B, Häusser H, Rüter A. Beschleunigt niedrig intensiver, gepulster Ultraschall die Heilung von Skaphoidfrakturen? *Handchirurgie, Mikrochirurgie, Plastische Chirurgie* 2000; 32:115-22

68. Muché JA, Efficacy of therapeutic ultrasound treatment of a meniscus tear in a severely disabled patient: a case report1 Arch Phys Med Rehabil 2003;84(10):1558–9
69. Nakamura T, Fujihara S, Yamamoto-Nagata K, Katsura T, Inubushi T, Tanaka E. Effects of low-intensity pulsed ultrasound on the expression and activity of hyaluronan synthase and hyaluronidase in IL-1<sub>α</sub>-stimulated synovial cells. Ann Biomed Eng. 2010;38(11):3363-70
70. Nakamura T, Fujihara S, Yamamoto-Nagata K, Katsura T, Inubushi T, Tanaka E. Low-intensity pulsed ultrasound reduces the inflammatory activity of synovitis. Ann Biomed Eng. 2011;39(12):2964-71
71. Ng CO, Ng GY, See EK, Leung MC. Therapeutic ultrasound improves strength of Achilles tendon repair in rats. Ultrasound Med Biol. 2003;29(10):1501-6
72. Nolte PA, Windt ESM De, Maas M, Dijk CN Van, Albers GHR. Het effect van ultrageluid op de botheling bij osteotomieën van de onderste extremiteit. Ned Tijdschr Traumat 1998; 6(5): 120-3
73. Nussbaum E. The influence of ultrasound on healing tissues. J Hand Ther. 1998;11:140-7
74. Nussbaum, EL (1997). "Ultrasound: to heat or not to heat - that is the question." Physical Therapy Reviews 2: 59-72
75. Oakley EM. Evidence for effectiveness of ultrasound treatment in physical medicine. Br. J. Cancer 1982;45Suppl.V:233
76. Ozgönel L, Aytakin E, Durmusoglu G. A double-blind trial of clinical effects of therapeutic ultrasound in knee osteoarthritis. Ultrasound Med Biol. 2009;35(1):44-9
77. Pilla AA, Mont MA, Nasser PR, Khan SA, Figueiredo M, Kaufman JJ, et al. Non-invasive low-intensity pulsed ultrasound accelerates bone healing in the rabbit. J Orthop Trauma 1990;4(3):246-53
78. Piravej K, Boonhong J. Effect of ultrasound thermotherapy in mild to moderate carpal tunnel syndrome. J Med Assoc Thai. 2004;87(Suppl 2):S100-6
79. Pope, GD, Mockett SP, et al. A survey of electrotherapeutic modalities: ownership and use in the NHS in England. Physiotherapy. 1995; 81(2): 82-91
80. Prentice, WE. Therapeutic Modalities in Sports medicine. 3rd edition, St.Louis, Mosby, 1994
81. Robertson VJ, Baker KG. A review of therapeutic ultrasound: effectiveness studies. Phys Ther. 2001;81(7):1339-50
82. Robertson VJ. Dosage and treatment response in randomised clinical trials of therapeutic ultrasound. Phys Ther Sport. 2002;3(3):124-33
83. Runday HK, Huizinga TWJ. Ultrageluid als behandeling voor het carpaletunnelsyndroom. Bed Tijdschr Geneesk. 1998;142(46):2534
84. Rubin C, Bolander M, Ryaby JP, Hadjiargyrou M. The use of low intensity ultrasound to accelerate the healing of fractures. J Bone Joint Surg Am 2001; 83-A(2): 259-70
85. Rutjes AW, Nüesch E, Sterchi R, Jüni P. Therapeutic ultrasound for osteoarthritis of the knee or hip. Cochrane Database Syst Rev. 2010;(1):CD003132
86. Rutten S, Nolte PA, Guit GL, Bouman DE, Albers GH, Use of low-intensity pulsed ultrasound for posttraumatic nonunions of the tibia: a review of patients treated in the Netherlands. J Trauma.

2007;62(4):902-8

87. Ryaby JP, Mathew J, Duarte-Alves P. Low intensity pulsed ultrasound affects adenylate cyclase and TGF- $\beta$  synthesis in osteoblastic cells. *Trans Orthop Res Soc* 1992; 17: 590
88. Sarrafzadeh J, Ahmadi A, Yassin M. The effects of pressure release, phonophoresis of hydrocortisone and ultrasound on upper trapezius latent myofascial trigger point. *Arch Phys Med Rehabil.* 2012;93(1):72-7
89. Shomoto K, Takatori K, Morishita S, Nagino K. Effects of Ultrasound Therapy on Calcificated Tendinitis of the Shoulder. *J Jpn Phys Ther Assoc* 2002; 5:7-11
90. Shu B, Yang Z, Li X, Zhang LQ. Effect of different intensity pulsed ultrasound on the restoration of rat skeletal muscle contusion. *Cell Biochem Biophys.* 2012;62(2):329-36
91. Smidt N, Assendelft WJ, Arola H, Malmivaara A, Greens S, Buchbinder R, van der Windt, DA, Bouter LM. Effectiveness of physiotherapy for lateral epicondylitis: a systematic review. *Ann Med.* 2003; 35(1):51-62
92. Sparrow, KJ, Finucane SD, et al. The effects of low-intensity ultrasound on medial collateral ligament healing in the rabbit model. *Am J Sports Med.* 2005;33(7):1048-56
93. Speed CA. Therapeutic Ultrasound in Soft Tissue Lesions. *Rheumatology* 2001; 40(12):1331-6
94. Takakura Y, Matsui N. Low-intensity pulsed ultrasound enhances early healing of medial collateral ligament injuries in rats. *J Ultrasound Med.* 2002; 21(3):283-8
95. Tascioglu F, Kuzgun S, Armagan O, Ogutler G. Short-term effectiveness of ultrasound therapy in knee osteoarthritis. *J Int Med Res.* 2010;38(4):1233-42
96. ter Haar, G. Therapeutic Ultrasound. *Eur J Ultrasound.* 1999; 9(1): 3-9
97. Tsai WC, Tang ST, Liang FC. Effect of therapeutic ultrasound on tendons. *Am J Phys Med Rehabil.* 2011;90(12):1068-73
98. Van der Windt DAWM, Van der Heijden GJMG, Van den Berg SGM, Ter Riet G, De Winter AF, Bouter LM. Therapeutic ultrasound for acute ankle sprains. *Cochrane Database Syst Rev.* 2006; 1
99. Walker NA, Denegar CR, Preische J. Low-intensity pulsed ultrasound and pulsed electromagnetic field in the treatment of tibial fractures: a systematic review. *J Athl Train.* 2007; 42(4):530-5
100. Ward AR, Robertson VJ. Comparison of heating of non-living soft tissue produced by 45 kHz and 1 MHz frequency ultrasound machines. *J Orthop Sports Phys Ther* 1996;23:258-66
101. Warden SJ, Bennell KL, McMeeken JM, Wark JD. Acceleration of fresh fracture repair using the Sonic Accelerated Fracture Healing System (SAFHS): a review. *Calcif Tissue Int* 2000;66:157-63
102. Warden SJ, Fuchs RK, Kessler CK, Avin KG, Cardinal RE, Stewart RL, Ultrasound produced by a conventional therapeutic ultrasound unit accelerates fracture repair. *Phys Ther.* 2006;86(8): 1118-27
103. Warden SJ. A new direction for ultrasound therapy in sports medicine. *Sports Med.* 2003;33(2): 95-107
104. Warden SJ, McMeeken JM. Ultrasound usage and dosage in sports physiotherapy. *Ultrasound Med Biol.* 2002;28(8):1075-80

105. Watson, T. The role of electrotherapy in contemporary physiotherapy practice. *Man Ther.* 2000; 5(3): 132-41
106. Wells PNT. Surgical applications of ultrasound. In: Nyborg WL, Ziskin MC. *Biological effects of ultrasound.* New York: Churchill Livingstone, 1985: 121-33
107. Wilkin LD, Merrick MA, et al. Influence of therapeutic ultrasound on skeletal muscle regeneration following blunt contusion. *Int J Sports Med.* 2004;25(1):73-7
108. Xavier CAM, Duarte LR. Treatment of nonunions by ultrasonic stimulation: first clinical applications. Meeting of the Latin-American Orthopedic Association, in conjunction with the 54th annual meeting of the American Academy of Orthopedic Surgeons; San Francisco; 1987 Jan 25
109. Yadav YK, Salgotra KR, Banerjee A. Role of Ultrasound Therapy in the Healing of Tibial Stress Fractures. *MJAFI,* 2008; 64(3): 234-6
110. Yang PF, Li D, Zhang SM, Wu Q, Tang J, Huang LK, Liu W, Xu XD, Chen SR. Efficacy of ultrasound in the treatment of osteoarthritis of the knee. *Orthop Surg.* 2011; 3(3):181-7
111. Yeung CK, Guo X, Ng YF. Pulsed ultrasound treatment accelerates the repair of Achilles tendon rupture in rats. *J Orthop Res.* 2006;24(2):193-201
112. Zeng D, Luo Q, Lin H, Zhang J, He C. The effect of therapeutic ultrasound to apoptosis of chondrocyte and caspase-3 and caspase-8 expression in rabbit surgery-induced model of knee osteoarthritis. *Rheumatol Int.* 2012; 32(12):3771-7
113. Zhou FH, Zhao HY. Acupuncture and ultrasound therapy for temporomandibular disorders. *Di Yi Jun Yi Da Xue Xue Bao.* 2004;24(6):720-1

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==== [www.sport-tec.de](http://www.sport-tec.de) =====

**Sport-Tec GmbH**  
**Physio & Fitness**  
Lemberger Str. 255  
D-66955 Pirmasens

**Tel.:** +49 (0) 6331 1480-0  
**Fax:** +49 (0) 6331 1480-220  
**E-Mail:** [info@sport-tec.de](mailto:info@sport-tec.de)  
**Web:** [www.sport-tec.de](http://www.sport-tec.de)

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Copyright: