THERA BEANS

ANLEITUNG

Art.-Nr.  29153 · Kategorie:  Therapie & Sinnestrainer

TheraBeans®

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Produktbeschreibung

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Manual instruction

Product-Description

TheraBeans® are an easy to use and highly efficient exercise-device for all the muscles of feet, lower leg, hands and forearms.

Developed from a training-method of high performance sports - the "corn-box" of Russian triple-jumpers, a sport with an extreme stress for the upper joint of the ankle and the stabilizing muscles of the foot - the TheraBeans® now offer a hygienically perfect utilization within the therapy of paretic muscles

Every muscle is trained with a simple set of exercises without being overloaded. At the same time a massage-effect with an increase of blood-circulation is achieved. Also strong sensorimotoric responses are created - the TheraBeans® therefore are an valuable tool for psychomotorical therapy too.

Priority of use is the rehabilitation of paretic muscles in upper and lower extremity (orthopedics, surgery, sport-rehabilitation, neurology) for a simple, efficient training of patients under own regimen.

Features

TheraBeans® are made of highly stable, PVC-free plastic. The used colors are cadmium-free and food-proof.

The two sizes are ideally tuned for the training and offer an even, well-dosed resistance to movements.

Resistance is generated by the mechanical friction of the beans when moved and is nearly neutral in direction - thus having no load-free phase - training in the beans is therefore always very easy in the beginning, though very exhausting for the muscles in a short time due to missing relaxation periods.

For a training the following amounts are needed:

Base-Set foot / forearm training
approx. 50 Liter / 30kg Case-size 45cm x 40cm x 40cm

Training both feet / both forearms
approx. 80 Liter / 50kg Case-size 60cm x 45cm x 40cm

The used case must be so large, that the TheraBeans and the body-part to be exercised can be taken into the volume even when moving.

TheraBeans® are easily cleaned and disinfected. Please take care for sufficient rinsing with clear water after disinfection, to avoid skin irritations when working later.

Regular Use

TheraBeans® are a training-device for a therapy of paretic muscles at various indications.

They are used in a suitable case, which is firm enough to cope with the weight of the beans and the load of the exercising person.

It must be large enough to provide ample space for the exercising person and permitting the exercise and the movement of the beans.

The case must be placed secured and slip-free. The exercising person must be able to maintain a comfortable and safe exercise-position.

Please take care for proper back-safe working when moving the filled case or use suitable transport-tools.

Preparations for the Exercises

Take off all jewelry and accessories (as possible), if needed cover irritated or wounded parts of skin.

Put arm or leg/foot, according to the desired exercise, slowly into the beans and work yourself downwards into the beans with small movements until your arm is covered with beans shortly above the elbow or your leg at least one hand below the knee.

Follow the adjoining exercise-instruction poster (from 15kg package on) or the correspondent instructions of your therapist for the specific exercises. Take care for minute execution of the instructions as there is a risk to be overloaded or wrongly loaded (problem e.g. at fractures).

Always start the exercises slowly and get accustomed to the movements and the resistance (warm up).

Never work into the pain, always remain within your personal range of tolerance - dangerous overload is thus avoided reliably.

Slowly increase the intensity of the exercise (velocity / range of motion) to your personal maximum possible strength and work out until you are no longer able to execute the movement properly or your muscles are exhausted.

Finish the training and take arm/leg out of the beans. Do not take up the work-out before finishing a relevant break (approx. 5 - 8 h).

Warnings

Children below the age of 3 are not allowed to work with and in the TheraBeans®, if not monitored.

TheraBeans® are not suitable to be eaten

TheraBeans® should only be used when properly instructed and on suggestion of a therapist

Do not use TheraBeans® at strong circulatory disorders (arterial occlusion etc.) and known thrombosis (risk of embolism).

Wounded or irritated skin parts must be covered and protected.

Exact execution of the exercise movements is essential at existing fractures or other instabilities in joints and bones to avoid dangerous overload. Please consult your therapist on this topic.

Technical Data

Material: Polypropylene
Density: 0.6 kg per liter

Cleaning

Suitable for standard detergents / disinfectants for regular cleaning (sometimes surface of beans will get dim after a time)

Conditions of storage, transport and use

Ideal conditions are temperatures above 0° Celsius. Avoid direct sun. The TheraBeans® contain no uv-stabilizer as this could only be achieved by recruiting plastic-softeners.

Temperature range is from 1°Celsius up to 50°Celsius, warmed up in oven or on heating (do not use microwave as the colour contains small metallic quantities)

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TheraBeans®

- Place the case with the TheraBeans® stable and slipfree and look for a suitable exercise position.
- According to the desired exercise: Sink your arm or leg so deeply into the TheraBeans® that the forearm is at least shortly above the elbow or the lower-leg one hand below the knee covered with beans and stay this deep within the whole exercise.
- Please cover and protect wounded or irritated parts of your skin.
- Children below the age of three must not work in or with the TheraBeans® without monitoring.

Exercises


- Always start the exercises slowly and get accustomed to the movement and the resistance (warm-up).
- Never work into the pain, always remain within your personal range of tolerance - dangerous overload is thus avoided reliably.
- Slowly increase the intensity of the exercises (velocity/range of motion) to your personal maximum possible strength and work out until you are no longer able to execute the movement properly of your muscles are exhausted.
- Finish the exercise and do not start again before having a relevant break - 5 to 8 hours normally sufficient.

Muscle-Exercises Foot - Extension/Flexion

movement:
Repeatedly lift up and press down your back of the foot.

time:
approx. 15 min.


Muscle-Exercises Foot - Pronation/Supination

movement:
Slightly lift up your foot and move in semicircle from outside towards inside. Lift up your foot in outward movement (pronation) and press down in inward motion (supination).


Muscle-Exercises Hand - Extension/Flexion

movement:
Repeatedly lift up your hand and press down again. Make a fist when pressing down, extend and spread your fingers when lifting up. Alternatively rotate your arm so that your hand stands sideways (the thumb looking upward)!


Muscle-Exercises Hand - Pronation/Supination

movement:
Rotate your hand and forearm out and inward (looking "into" the fist = supination, looking "onto" the back of the hand = pronation). Extend and spread your fingers in inward rotation, make a fist at outward rotation.

Muscle-Exercises Foot - Rotation

movement:
Lift up your foot and make a complete circular movement. Extend your foot in the inward rotation (supination), in outward rotation (pronation) lift your foot.


Sensomotor Exercises (without picture)

- Grab just one particular bean (with your hand or foot)
- You determine a certain number of beans (up to 40 pieces). Try to grab these at once and control the result.
- Differentiate between the two sizes of beans. Then grab alternately just the big or small ones only or follow instructions to grab certain compositions of large and small beans.
- Hide something in the pile of beans and try to find the object (the smaller the object, the more difficult this task is).
- Let somebody hide an unknown object in the beans, try to trace this again and recognize the object by touching.
- Move the object forwards in the beans by pushing it smoothly - do not grab it!

Muscle-Exercises Foot - Combined Motion

movement:
Do not lift your foot! "Walk" sideways in the beans from left to the right and back. Only move your foot, do not push with the lower leg (dangerous overload)!

Muscle-Exercises Hand - Combined-Motion

movement: (without picture)
Move your hand in the beans, as if you stir something. Repeatedly change direction.

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copyright Jakobs GmbH
Sport-Tec GmbH