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## A healthy foot -

### the fundamental base of our body

- Carries our body weight
- Regulates our equilibrium
- Provides a variety of movement possibilities

With healthy feet we can walk, jump, run, dance and hop. To fulfil all these duties, the foot depends on an elastic arch structure.

### Foot problems

About forty percent of the population suffers from foot problems of which about two thirds are caused by movement habits and wrongly applied pressure. This results in skew-, flat- or splayfeet. These defective positions are often the beginning of a concatenation of bigger foot problems and the whole leg axes gets out of plumb. This means inevitably strain of knee joints and hip joints, as well as of backbone.


### Foot screwing (anatomy of foot)

Our foot is a sensitive entity of 26 bones, 27 joints, 107 ligaments and 32 muscles and tendons. Dr. Christian Larsen discovered that nature has built our feet after the spiral principle, which means: the fore foot turns inwards (prone extension), the back foot (heel) outwards (supine flexion). The so originating length- and crosswise vault is supported by three sphenoid bones. Wrong foot position causes loss of foot curvature and the foot collapses in itself. Consequently toe and metatarsal bones will be displaced.

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
## A healthy foot -

### the fundamental base of our body

In order for the foot screwing to dynamically push off and absorb the shock, the foot must be adequately loaded. Dr. Larsen developed a training and movement specifically for this, called spiral dynamic.

**Pro-Pedes**, the foot-leg-axe equipment of **Pedalo**, offers many exercises to create and maintain a healthy foot.

We wish you a lot of fun and success with the **Pro-Pedes**

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Antirutsch-Matte Art.-Nr. 628099

## Pro-Pedes

This foot-leg-axle equipment is ideal to learn and train the foot-screwing for a healthy foot. While standing in an accurate heel position and having big-toe-joint/ball contact the longitudinally arch is trained and at the same time a fine and straight leg-axle is stabilized. Used with the flat side on top Pro-Pedes is an effective therapy roundabout to improve posture and movement coordination, characterized by best rolling. Proven for children, adults, senior citizens as well as after injuries in postural system.

Item no. 143 160

### Specification:

Solid beech, arched on one side, approx. 270x49 mm, linen bag and exercise instructions included.

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## A) Activation of foot screwing (FS)...

✓ Do this exercise over the convex part of the Pro-Pedes to support the longitudinally arch (counteracts skew and flat feet).

### FS 1 while sitting:

**Position:** Sitting on a stool, put both feet next to each other. Make sure both knees and both feet point forward. Place your left foot concentric on the Pro-Pedes.

**Practice:** Lift up your left foot, lower it slowly and place it concentric on the Pro-Pedes. Make sure your heel moves to the back and down and the big toe joint moves forward/down. During the next 30 seconds intensify the contact in the external third of the heel at the same time as in the big toe joint.

**Duration:** Repeat 10 times.

### FS 2 while standing:

**Position:** Same as FS1 but in standing position. Please be careful, practice regularly but only briefly until your foot has accustomed to the foot screwing.

**Practice:** Lift up your left foot, lower it slowly and place it concentric on the Pro-Pedes. Make sure your heel moves to the back and down and the big toe joint moves forward/down. During the next 30 seconds intensify the contact in the external third of the heel at the same

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multiple posture and movement coordination

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## A) Activation of foot screwing (FS)...

### FS 3 in one-leg-stand:

**Position:** Place your foot concentric on the Pro-Pedes, heel moves to the back down, big toe joint moves forward down.

**Practice:** Keep standing on one leg and balance while standing with straight leg and leg-axle.

Make sure toes and kneecap point forward.

**Duration:** Balance up to 2 minutes, repeat 3 times.

### FS 4 in one-leg-stand, free leg is moving:

**Position:** Place your foot concentric on the Pro-Pedes, heel moves to the back down, big toe joint moves forward down.

**Practice:** Do with your free leg:

- draw an 8 in the air
- write your name in the air
- throw a ball to against a wall or to a partner

**Duration:** Balance up to 2 minutes, repeat 3 times.

### FS 5 while walking:

**Position:** Standing in front of the Pro-Pedes.

**Practice:** Walk over the Pro-Pedes with your heel coming down on the outer third of it, the mid of your foot on the top of it and the big toe joint/ball landing secure on the Pro-Pedes. Make sure your kneecap keeps pointing forward.

**Duration:** Walk 10 times over the Pro-Pedes while paying attention to a straight leg-axle.

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## B) Activation of foot lateral arch on forefoot (FLA)...

✓ This exercise is done over the convex part of the Pro-Pedes to activate the lateral arch (counteracts the splayfeet).

### FLA 1 while sitting:

**Position:** Sitting on a stool, put both feet next to each other. Make sure both knees and both feet point forward. Place your left foot on the highest point of the Pro-Pedes. Your heel stands on the edge.

**Practice:** Form your lateral arch over the convexity of the Pro-Pedes, by moving your knuckles of small and big toes inwards and outwards. Hereby toes stay relaxed, your heel stays down.

**Duration:** Repeat 10 times.

### FLA 2 while standing:

**Position:** In standing position: Put both feet next to each other. Make sure both knees and both feet point forward. Place your left foot on the highest point of the Pro-Pedes. Your heel stands on the edge.

**Practice:** Form your lateral arch over the convexity of the Pro-Pedes, by moving your knuckles of small and big toes inward and outward. Hereby toes stay relaxed, your heel stays down.

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## B) Activation of foot lateral arch on forefoot (FLA)...

### FLA 3 while standing in one-leg-stand:

**Position:** In one-leg-stand. Place your foot with lateral on the highest point of the Pro-Pedes. Your heel stands on the edge.

**Practice:** Keep standing on one leg and balance while standing with straight leg and leg-axle. Make sure toes and kneecap keep pointing forward.

**Duration:** Balance up to 2 minutes, repeat 3 times.

### FLA 4 while standing in one-leg-stand, free leg is moving:

**Position:** In one-leg-stand. Place your foot with lateral on the highest point of the Pro-Pedes. Your heel stands on the edge.

**Practice:** Do with the free leg:

- swing forwards and backwards
- draw an 8 in the air
- write your name in the air
- throw a ball against a wall or to a partner

**Duration:** Balance up to 2 minutes, repeat 3 times.

### FLA 5 while walking:

**Position:** Stand in front of the Pro-Pedes.

**Practice:** Walk with your heel on the outer edge and with the lateral on the highest point of the Pro-Pedes.

**Duration:** Walk 10 times over the Pro-Pedes while paying attention to a straight leg-axle.

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## C) Stretching/lengthening of the Achilles' tendon (AT)...

✓ Improves mobility in your ankle which enables you to land on your heel better when walking. This leads to an upright body posture.

### AT 1 while standing:

**Position:** Place the Pro-Pedes approx. 40-50 cm from the wall onto an Anti-slip-mat. Position your forefoot on the highest point of the Pro-Pedes. Toes and kneecap point forward. Support yourself with both hands in front of you onto the wall.

**Practice:** Push yourself up onto your toes as easy as possible. Afterwards - in slow motion - lower your heel deeply. Remain with your heel down for 30 seconds, then push yourself up - in slow

motion - onto your toes again to then repeat the process of lowering your heel deeply, remaining for 30 seconds and so on ...

**Duration:** Repeat 10 times

### AT 2 while standing in one-leg-stand:

**Position:** Same as AT1

**Practice:** Push yourself up onto your toes as easy as possible. Afterwards - in slow motion - lower your heel deeply. Remain with your heel down for 30 seconds, then push yourself up - in slow motion - onto your toes again to then repeat the process of lowering your heel deeply, remaining for 30 seconds and so on ...

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## D) Leg-Axle stabilisation (LA)

Strengthening of the leg musculature -  
by balance training ...

### LA 1 while standing:

**Position:** Stand with both feet on the flat side of the Pro-Pedes. Toes and kneecaps point straight forward, knees are slightly bent.

**Practice:** Balance so that the Pro-Pedes always remain in a horizontal position.

**Duration:** Balance up to 2 minutes, repeat 3 times.

### LA 2 while standing in one-leg-stand:

**Position:** Stand with one leg on the flat side of the Pro-Pedes. Toes and kneecap point straight forwards, knee is slightly bent.

**Practice:** Balance so that the Pro-Pedes always remain in a horizontal position.

**Duration:** Balance up to 2 minutes, repeat 3 times.

### LA 3 while standing in one-leg-stand, free leg is moving:

**Position:** Stand with one leg on the flat side of the Pro-Pedes. Toes and kneecap point straight forward, knee is slightly bent.

**Practice:** With the free leg:

- swing forwards and backwards
- draw an 8 in the air
- write your name in the air
- throw a ball against a wall or to a partner.

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- gehen
- steigen Treppe

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## E) Posture exercises (PE)...

✓ To develop a feeling for the correct upright posture.

### PE 1 while standing:

**Position:** Standing on both legs you put the Pro-Pedes with the flat side on your head.

**Practice:** Balance the Pro-Pedes on your head.

**Duration:** Up to 2 minutes, repeat 3 times.

### PE 2 while standing in one-leg-stand:

**Position:** Standing on both legs you put the Pro-Pedes with the flat side on your head.

**Practice:** Balance the Pro-Pedes on your head.

**Duration:** Up to 2 minutes, repeat 3 times.

### PE 3 while standing in one-leg-stand, free leg is moving:

**Position:** Standing on one leg, put the Pro-Pedes with the flat side on your head.

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**Duration:** Balance up to 2 minutes.

### PE 4 while walking:

**Position:** Place the Pro-Pedes with the flat side on your head.

**Practice:** Balance the Pro-Pedes on your head and

- walk
- go upstairs and downstairs or imitate this type of movement.

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## Further exercises

✓ **Used with the flat side up Pro-Pedes is an effective therapy roundabout to improve body posture and movement coordination, with great rolling motion. It has proven successful for children, adults, senior citizens as well as after injuries in postural system.**

**Position:** Standing with both legs on the Pro-Pedes, knees are slightly bent, looking ahead.

**Practice 1:** Try to keep the platform of the Pro-Pedes as horizontal as possible. Duration: 30 sec.

**Practice 2:** based on practice 1, turn your head alternating left/right for 5 seconds each. Repeat 4 times. Then put your head back (look up) for 5 seconds, afterwards look ahead for 5 seconds. Repeat 4 times.

**Practice 3:** Squat slowly deeply and keep this position for 2 seconds, straighten up slowly and keep this position for 2 seconds. Duration: Repeat 4 times.

**Practice 4:** Wiggle with the Pro-Pedes forward/backward, left/right.

**Position:** Standing in one-leg-stand, knee is slightly bent.

**Practice 1:** Try to keep the platform of the Pro-Pedes as horizontal as possible. Duration: 30 sec.

**Practice 2:** Swing your free leg forward/backward.

**Duration/practice:** 30 seconds each, then change leg.

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## Important information

Please read and follow!

The Pedalo-Pro-Pedes is made of untreated massive wood. The wood adapts itself to the surrounding climate constantly and thereby changes his size (wood lives). If this process happens too fast it results in cracking. Hence make sure that Pro-Pedes will not placed on a heating or in the direct sun at any time. Also permanent storage on a ground with under-floor heating can lead to the cracking.

Please store the Pro-Pedes in the linen bag at an evenly air-conditioned place after using, e.g. in a wardrobe.

In case you wipe the Pro-Pedes moistly sometime, please rub it dry afterwards immediately with a cloth and let it dry at an evenly air-conditioned place.

Dear customer,

The Pedalo-article you have chosen is made of natural wood. Differences in colour, structure, grain, small included knots or scars prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

**Surface:**

Delivered untreated. If you wish you can oil the beechwood disc with every natural oil or just leave it raw – as is. To have the disc slightly gripped you might rub the surface with a little water one time, so the wooden fibers straighten up and the surface will have more grip.

**Please note:**

For the use of this device a certain degree in motor abilities and circumspection is assumed. In case of doubt first use device only with assistance.

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## Further products for the „healthy foot“

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
### Pedalo® Massage-mat:

3 minutes each day are enough not only to release your foot musculature and consequently the deep musculature of the whole body, but also to strengthen it. Your feet become injured, the forced up blood circulation activates your metabolism and hence contributes to a faster regeneration in sports and every day life.

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- ✓ Strengthening of foot and lower leg musculature

### Pedalo® PhysioFlip:

PhysioFlip is mainly used for the improvement of the effectiveness of different muscle groups, e.g.:

- lower leg musculature (knee-joint rotators)
- thigh musculature (hip-joint rotators)
- upper arm and shoulder musculature (shoulder-joint rotators)
- with foot complaints: flat feet, hollow feet, deformations of the foot, overpronation

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
### Pedalo® Footboards moveable fulcrum:

The Footboards moveable fulcrum can be used individually and allow varied exercises for foot-leg-axle training to improve coordination and stabilisation, strengthening and stretching of the foot to erect the horizontally foot arch. Your stability and walking safety will be improved noticeably.

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- ✓ Sensitization of reflex zones
- ✓ Relaxation of deep musculature

### Pedalo® Pedasan:

Exercising barefoot on Pedalo-Pedasan helps to strengthen the foot musculature and to boost the neglected foot vault. The reflex zones of the foot soles are sensitized what activates the reaction of muscles, tendons and nerves. Exercises with Pedalo-Pedasan additionally train balance, reaction and coordination.

**Sport-Tec**  
Physio & Fitness

Lemberger Straße 255  
D-66955 Pirmasens

Tel.: 06331/1480-0  
Fax: 06331/1480-220

info@sport-tec.de  
www.sport-tec.de