BLACKROLL®

NEEDLEROLLER

Instruction Manual

EN | DE | FR | NL | IT | PL | CN

Product description

Congratulations on your new BLACKROLL® NEEDLEROLLER! Stimulating the superficial fasciae with the NEEDLEROLLER has a variety of positive effects. The nerve endings in the skin and the upper fascia layers are stimulated, which has a circulation-promoting effect and can reduce sensitivity to pain. The NEEDLEROLLER is particularly suitable for stimulating reflex and acupuncture zones. When used correctly, it can regulate the production of collagen and elastin and optimally prepare the tissue for subsequent therapeutic treatments.

Safety instructions

- + Avoid excessive pressure: If the stimulation is too intense or redness and irritation occur, reduce the pressure or take a break.
- + Do not use if you are allergic to nickel: The NEEDLEROLLER contains metals, including nickel. It should not be used if you are allergic to nickel.
- + Consult a doctor before use if you suffer from acute or chronic illnesses, such as cardiovascular problems, acute infectious diseases, varicose veins or skin diseases.
- + Do not use on open wounds or fresh injuries.
- + Storage: Keep the NEEDLEROLLER in a dry place and avoid direct sunlight and high humidity (over 80%).

Contraindications

Do not use the BLACKROLL® NEEDLEROLLER if you suffer from acute infectious diseases, cardiovascular diseases, thromboses, varicose veins or other serious health complaints without first consulting your doctor. The roller should not be used by people with a nickel allergy.

How to use

1 Preparation

- + Gently roll the BLACKROLL® NEEDLEROLLER over the skin, adjusting the pressure to suit your individual needs. The pressure should be comfortable and not cause any pain.
- + Use the NEEDLEROLLER to stimulate the free nerve endings in the skin and to achieve a circulation-boosting effect.

2 Use

- + Use the NEEDLEROLLER to stimulate reflex and acupuncture zones.
- + Start with light pressure and increase as needed. Roll back and forth to stimulate the desired zones.
- + You can place a thin cloth between the NEEDLEROLLER and the skin to reduce the intensity of the stimulation if it is too strong.

3 After use

- + Clean the NEEDLEROLLER after each use with soap and water and dry it thoroughly with a hair dryer to prevent rusting.
- + Store the NEEDLEROLLER in a dry and well-ventilated place.

Technical data

Dimensions:	26 cm x 9 cm x 6.5 cm
Weight:	180 g
Needles:	steel (25%), steel with zinc coating (12.5%), nickel (25%), copper TM (25%), copper with silver coating (12.5%)
Handle:	ash wood
Base and plug:	rubber
Frame:	steel

Care instructions

- + Clean the BLACKROLL® NEEDLEROLLER regularly with soap and water and a soft brush.
- + Dry the roller thoroughly with a hair dryer after cleaning to prevent rusting.
- + Store the NEEDLEROLLER at a temperature between 10 and 40 degrees Celsius and a humidity of less than 80%.
- + Avoid storing it in direct sunlight or in humid environments.

Possible applications

- + Pain relief: The NEEDLEROLLER can help to interrupt pain cycles and reduce sensitivity to pain.
- + Stimulation of the skin: Use the roller to stimulate the superficial fascia layers and promote blood circulation.
- + Preparation of the tissue: The NEEDLEROLLER prepares the skin and tissue for subsequent therapeutic treatments.

Scope of delivery

1x BLACKROLL® NEEDLEROLLER
1x BLACKROLL® NEEDLEROLLER Manual

Warranty and liability

The BLACKROLL® NEEDLEROLLER comes with a two-year warranty. Damage caused by improper use or improper storage (e.g. rust formation due to moisture) is excluded from the warranty.

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BLACKROLL AG
WSVK Oederan GmbH