www.sport-tec.de

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**Propriomed, 100 cm**

**Art.Nr.: 23475**

**▶▶▶ [zum Produkt ...](#)**

**Propriomed**

**▶▶▶ [zur Kategorie...](#)**

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
# PROPRIOMED

Exercise instructions for postural training according to Dr. Rašev (concept by C. Puta and M. Herbsleb)

English

# PROPRIOMED: How does it work?

## General description of how it works

The physical action of the PROPRIOMED consists in the fact that the movements of the person exercising and the specific oscillating properties of the pole create cyclically triggered oscillations with defined frequencies, which are carefully dampened. The neuro-physiological effect comes to light when efferent impulses, with dosed stimulus density, are targeted at the muscles, which in turn leads to a stimulation of synergistic muscle activation (coactivation). This effect is crucial with respect to any stability requirements of the motor system when aligning a vertical body against gravity.

It is through dampening and individually adjusted regulation of the pole's own natural frequency as well as through an amplitude that is set independently that both the level of sensitivity and control performance of the motor system can be optimized in a targeted manner.

## Vibrational properties of the PROPRIOMED

The PROPRIOMED's action can be influenced manually on the basis of two physical properties.

These are:

- Changes in its natural frequency by manually adjusting the frequency controllers
- Changes in the amplitude by employing varying degrees of power transmission

## Frequency

The natural frequency with which the PROPRIOMED pole can be stimulated to vibrate can be varied by means of manually adjustable controllers between 2.5 and 7.5 Hz, depending on the device. This allows for a seamless dose that can be adjusted for the individual. These adjustment options provide the following benefits in terms of preventive training and therapy:

- Global and regional activation of the neuro-muscular system: individual regions (thoracic girdle) are particularly stimulated when the whole body (global) is stabilized
- Individual, seamless adjustment in line with the performance level as well as the training and therapeutic objective
- Asymmetrical adjustment of the frequency controllers



Essentially, the following needs to be taken into consideration:

"The shorter the pole length, the more difficult the exercise."

In other words, to make the exercise as difficult as possible, the frequency controllers would be right inside the grip.

Note that the number in the device designation corresponds to the length of the device (in centimetres).

An exception of this principle is the PROPRIOMED 100. Special therapeutic requirements were instrumental in the development of this "lightweight". The PROPRIOMED 100 is ideal, due to its low own weight and frequency band, for treating people right after shoulder surgery, reconstruction of ligaments in the shoulder joint and acromio-clavicular joint (joint between the acromion – at shoulder level, outer end of the slightly sloping, sturdy crest at the back of the shoulder blade – and the collar bone). The PROPRIOMED 100 is also suitable for children and adolescents due to its low own weight.

Manual adjustment of frequency controllers

The following graphic illustrates the natural frequency of the PROPRIOMED in relation to the length of the device. This allows for the setting of six levels of difficulty.


Table showing the setting of levels of difficulty


### Amplitude

The vibration's amplitude, which is also continuously adjustable, supports the device's multivariable use in preventive training, fitness and therapy.

At a maximum, the device-dependent oscillation amplitude should not be exceeded. Any increases in intensity should be implemented by adjusting the frequency controllers or changing to the next smaller device.

	Length in cm	Weight in g	optimal amplitude	max. amplitude
□	100	504	30 cm	40 cm
■	190	1081	60 cm	70 cm
■	170	1035	50 cm	60 cm
■	150	994	40 cm	50 cm
■	130	947	30 cm	40 cm

**Do not exceed the maximum amplitude! Exceeding the amplitude will reduce the useful life of the device!**




## Oscillating planes 1 and 2

The PROPRIOMED can be used on two different oscillating planes - oscillating plane 1 and oscillating plane 2.

### Oscillating plane 1


With respect to this plane, the PROPRIOMED is stimulated to oscillate as an extension of the lower arm (i.e., in the direction in which the lower arm points). Movement is "forward and back". 

### Oscillating plane 2

As for oscillating plane 2, the PROPRIOMED is moved against the extension of the lower arm. This corresponds to a change from oscillating plane 1 by 90°.

Depending on the respective hand position, movement is "left and right"  or "up and down". 

## Information regarding the exercise

Below, you see an example to explain an exercise, including the legend. The symbols shown here are used in all exercise sheets and the poster (page 21).

**S** = Start

**E** = End

= oscillating plane 1: the PROPRIOMED is moved "forward and back".

= oscillating plane 2: the PROPRIOMED is moved "left and right" or "up and down".

Legend

Please read the following table carefully. You should internalize the terms used here and their explanations to ensure that exercises are carried out successfully.

Term	Explanation
<b>Exercise Unit (EU)</b> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EU 1 5 - 15 sec</div>	An exercise unit represents uninterrupted exercise time. This should be between 5 to 15 seconds. A good yardstick is counting "in your mind" to 10.
<b>Repeated Exercises (RE)</b> <div style="display: flex; align-items: center; gap: 10px;"> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EU 1</div> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EU 2</div> <div style="font-size: 20px;">...</div> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EU 10</div> </div>	A repeated exercise is equivalent to the one-time repetition of the exercise unit. The maximum number of repetitions should not exceed 5 to 10.
<b>Exercise Cycle (EC)</b> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EC 1 20 - 30 min</div>	An exercise cycle corresponds to the one-time completion of all exercises including 5 to 10 exercise repetitions each. The maximum number of exercise cycles should not exceed 2. An exercise cycle corresponds to an overall exercise time of approx. 20 to 30 minutes.
<b>Exercise Break (EB)</b> <div style="border: 1px solid black; padding: 2px; display: inline-block;">EB 5 - 10 sec</div>	Each exercise unit should be followed by an active exercise break of 5 to 10 seconds by shaking out and loosening arms and hands.

## Basic structure of an exercise lesson


Your exercise lesson should be structured according to the principle illustrated in the diagram below. From a therapeutic point of view, this structure applies only in part. Specific recommendations can be found in the therapy program.

### Warming-up

Ideally, you should do general and specific warming-up exercises prior to working out with the PROPRIOMED. This should take you 5 to 10 minutes.

General warming-up consists of specific exercises that involve the whole body. These include, for example, bicycle riding, jump roping, walking or physical exercises.

Specific warming-up is about mobilizing the spine and the large extremity joints, in a targeted manner, such as hunching your back, circular movements with your arms, legs and feet (preferably from behind outward) as well as pendular movements.

### Criteria for stopping exercise and indications of excessive strain

Please observe the following criteria for stopping the exercise and indications of excessive strain:

- Pain
- Synkinesis (involuntary movement) of hands, thoracic and/or pelvic girdle
- The person can no longer keep the PROPRIOMED oscillating
- Cramping up of muscles
- Maximum exercise time (15 sec.) has been reached

To ensure the effectiveness and specificity of each exercise, it is especially important to observe and assume the best possible posture. On the one hand, this helps to relieve the passive structures, while, on the other hand, it optimizes the strain on the muscles targeted under this program.

Ideal posture is characterized by the deliberate segmental alignment of one's head, thoracic girdle, spine, pelvis and feet while standing up against gravity.

## General optimal posture

Please allow yourself sufficient time to implement the following steps towards optimizing your general posture consciously.

### Stand (general)

- Shoulder-wide stand (both feet)
- Hip joint and knee joint slightly bent (approx. 10°)
- Feet turned outward, approx. 10-15°
- Back of your head, thoracic girdle and buttocks form a straight line

Posture

## Segmental alignment

### Head position

- look forward and slightly downward
- double chin (move back of the head backwards, without lifting the chin)

### Upper body

- raise chest (raise sternum slightly forward and up ["showing off one's medals"])
- thoracic girdle rests on the chest
- let shoulder blades drop in the back (like putting your hands in your pockets)
- relax abdominal wall (primarily abdominal breathing)

### Leg axis

- slightly bend hip and knee joints (bend erect trunk slightly forward)
- legs slightly rotated outward
- note: thigh, calf and foot form one plane (no O or X-legged position)

### Feet

- balanced foot load (load distribution between forefoot and backfoot = 1/2 : 1/2)
- barefoot, if possible
- "short foot" according to JANDA: stretch toes, then splay them out and press them down on the floor, so that the small muscles along the foot's interior side become relaxed

## Notes

Our everyday lives are dominated by lack or monotony of movement (long drives in the car, working on the computer, etc.). Activities performed while standing up or sitting down that involve recurring, monotonous and delicate static or dynamic movements (mostly forward and towards the inside of the body) make up the bulk.

We often tend to assume a "hunched-over" posture, which leads to the development of certain negative factors and increases in misplaced strain on the skeletal system. Functional impairment of muscles and vegetative functions (e.g., breathing) can also be observed.

By assuming the posture described in this brochure and complying with the general and specific information during the exercises, it is possible to counteract this monotonous, everyday strain.

Depending on the style of exercise, the way in which you hold the device and the positioning of the PROPRIOMED, the most important groups of muscles are trained and exercised through coordinative interaction.

*Pain affecting the motor system is often rooted in functional disorders of the neuro-muscular system and is less often the result of structural failure!*

### General information regarding the use of the PROPRIOMED:

Following this information will help you learn the use of the PROPRIOMED systematically.

- the impulse is transferred from the body to the PROPRIOMED; it is triggered by moving your hand rhythmically forward and backward and/or up and down
- keeping the PROPRIOMED oscillating by repeated short impulses
- the oscillating amplitude (the swing of the PROPRIOMED) varies depending on the available potential force and the type of exercise. Do not exceed the maximum amplitudes of the PROPRIOMED under any circumstances (see p. 05)
- movement amplitudes in the thoracic and pelvic girdles should be kept to a minimum
- following each exercise, shake out and loosen your hands and arms (loosening of finger flexors)
- recommended exercise repetitions and breaks should be observed


The ability to make the PROPRIOMED oscillate and keep it oscillating is the first and, initially, the most important challenge when using this device!

At the same time, make sure that your hand movement is kept to a minimum; this way, you will have to maximize the power and energy from your body centre in order to achieve the desired effect for your muscles.


Please see the detailed exercise description for exercise-specific information.


Exercise 1


Exercise 2

Exercise 6

Exercise 3

Exercise 5


Exercise 4

Program

# Exercise 1

Plane 1

## Posture

For general information on optimal posture, please see p. 09.

Focus on the following steps. Trace them in your mind at first and then put them into practice.

## Type of grip and positioning of the PROPRIOMED

- two-hand grip/hands crossed
- PROPRIOMED is held vertically in the centre in front of your body

## Execution

- impulse transfer by moving your hand forward and backward
- while maintaining the oscillation, turn your upper body left by 5-10°, back to the centre and then to the right
- during the exercise, keep your pelvis steady
- let your head follow the rotation of the upper body

## Special notes


- oscillation amplitude only 20 to 30 cm
- change hand position after each EU (upper hand down; lower hand up)


## Strain parameters

EU: 5-10

Break: 5-10 sec.


## Targeted body regions

## Exercise 2

Plane 1

### Posture

Modify the general optimal posture as follows:

- knee joint: bent approx. 110°
- upper body: bent forward approx. 30°

Focus on the following steps. Trace them in your mind at first and then put them into practice.

### Type of grip and positioning of the PROPRIOMED

- two-hand grip / hands side by side (thumbs are pointed at each other)
- PROPRIOMED is held in parallel in the centre in front of your body, slightly above your knee joints

### Execution

- impulse transfer by moving your hand forward and backward
- during the EU, the starting position remains unchanged
- keep your pelvis steady throughout the entire exercise

### Special notes

#### Leg axis

- kneecaps do not "cover the tips of your feet"
- positioning of the knees: visually, an extension of the 2nd toe

#### Feet


- load distribution regarding the foot (forefoot to backfoot = 2/3 : 1/3)

### Strain parameters

EU: 5-10




Break: 5-10 sec.

### Targeted body regions

## Exercise 3

Plane 1

### Posture

For general information on optimal posture, please see p. 09.

Focus on the following steps. Trace them in your mind at first and then put them into practice.

### Type of grip and positioning of the PROPRIOMED

- two-hand grip/hands side by side (thumbs are pointed at each other)
- PROPRIOMED is held in parallel in the centre in front of your body, at the level of your navel
- both arms are bent at the elbows (approx. 20-30°)

### Execution

- impulse transfer by moving your hand forward and backward 
- while maintaining the oscillation, lift the PROPRIOMED over your head in a semi-circular motion (elbows approximately reach up to the ears)
- then lower it to the level of your navel
- this sequence corresponds to one EU

### Special notes


- Do not pull shoulders up during the exercise and keep your arms bent

### Strain parameters


EU: 5-10

Break: 5-10 sec.

### Targeted body regions


## Exercise 4

Plane 1

### Posture

For general information on optimal posture, please see p. 09.

Focus on the following steps. Trace them in your mind at first and then put them into practice.

### Type of grip and positioning of the PROPRIO-MED

- single-hand grip
- exercise-side hand grips the PROPRIO-MED from above at the level of the navel
- back of the hand points up and outward
- PROPRIO-MED is held diagonally in the centre in front of your body
- elbow is slightly bent
- free palm is pointed forward, with the fingers stretched and splayed out

### Execution

- impulse transfer by moving your hand forward and backward 
- while maintaining the oscillation, move the PROPRIO-MED outward-up-back in semi-circular motion
- in the end position, the back of the hand is pointed down
- hand and shoulder level form one line
- always keep the elbow below that notional line

### Special notes


- change exercise hand after each EU

### Strain parameters

EU: 5-(10)



Break: 10 sec.

### Targeted body regions

## Exercise 5

Plane 1

### Posture

Modify the general optimal posture as follows:

#### Free leg

- hip joint bent, approx. 30°; slight outward rotation of the thigh
- calf vertically down (parallel to support leg)
- lift toes and forefoot; imagine climbing stairs

Focus on the following steps. Trace them in your mind at first and then put them into practice.

### Type of grip and positioning of the PROPRIO-MED

- two-hand grip/hands side by side (thumbs are pointed at each other)
- PROPRIO-MED is held parallel in the centre in front of your body, at the level of the navel

### Execution

- impulse transfer by moving your hand forward and backward 
- the positions of the arms, the free leg and support leg remain unchanged throughout the entire exercise

### Special notes

- change support leg and free leg after each EU

#### Abort the exercise if:

- you need to stabilize your one-legged stand by jumping
- the knee of the support leg tilts inward
- the inside arch of the support leg tilts inward


Once the exercise has been aborted, continue with the other side.



### Strain parameters

EU: 5-(10)

Break: 10 sec.

### Targeted body regions

### Posture

For general information on optimal posture, please see p. 09.

Focus on the following steps. Trace them in your mind at first and then put them into practice.

### Type of grip and positioning of the PROPRIOMED

- two-hand grip/hands side by side (thumbs are pointed at each other)
- PROPRIOMED is held parallel in the centre in front of your body, immediately below the sternum
- both arms are bent at the elbows, approx. 20-30°

### Execution

- impulse transfer by moving your hand up and down 
- hands and arms do not change positions in space
- hands and elbows remain almost in a single line
- maintaining the oscillation in this position

### Special notes


- during the exercise, do not pull up your shoulders and maintain the bent position of your arms
- do not bend your wrists
- maintain segmental alignment during the exercise (e.g., your upper body must not sway forward or back)

### Strain parameters

EU: 5-10

Break: 5-10 sec.

### Targeted body regions

### **Short foot according to JANDA**

The short foot according to JANDA is the basic element in the concept of sensorimotor facilitation according to JANDA-VARVROVA, which is based on the physiotherapy according to Mensen-Dieck. The objective of this exercise is to systematically activate the foot's inside muscles. This is to increase the afferent flow to the central nervous system.

### **Neuro-muscular**

Neuro-muscular refers to the continuous interaction between muscles and nerves. "Movement is merely one expression of the control work performed by the CNS!" (Henner, 1946)

### **Surface electromyography**

For the purposes of surface electromyography, electrodes are attached to the skin to measure and evaluate the extent of muscular activation by way of potential recording.

### **Segmental alignment**

Segmental alignment is defined as the positioning of these body segments: thoracic girdle, pelvic girdle, head and spine in relation to each other. This segmental alignment is achieved by intersegmental muscles.

In the event of a postural segmental instability, the polysegmental muscles take over the function of the intersegmental muscles, which frequently results in pain due to the excessive strain.

### **Intersegmental muscles**

(Inter-)segmental muscles cover only a joint or segment of the spine and are primarily responsible for static work.

### **Polysegmental muscles**

Polysegmental muscles exist in two layers and are of varying length. Shorter muscles cover 4-6 segments of the spine, while the longer muscles are superficial and cover more than 6 segments of the spine. These muscles are primarily responsible for motor work.

### **Synergistic muscle activation**

See "Synergistic coactivation"

### **Synergistic coactivation**

Synergistic coactivation means that, in a simplified joint model, muscle pairs with opposite functions (flexors and extensors, external rotators – internal rotators, abductors – adductors) are activated simultaneously. This coactivation ("co-" from the Latin word meaning "with") is of the utmost importance to the performance of stabilizing (static) tasks. The intensity of such activation varies continuously depending on the task at hand.

## Information regarding exercise poster

All the exercises as well as the information on how to do the exercises and on posture are summarized again in the exercise poster attached to help you with the practical use of the PROPRIOMED.

Before you rely exclusively on the poster for your exercises, make sure you spend enough time first to work through the exercises in your mind and practically.

**Sport-Tec**  
Physio & Fitness

Lemberger Straße 255  
D-66955 Pirmasens

Tel.: 06331/1480-0  
Fax: 06331/1480-220

info@sport-tec.de  
www.sport-tec.de