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Paraffinbad Para Pro
inkl. 2,5 kg Wachs, 3 Liter

Art.Nr.: 29142

▶▶▶ [zum Produkt ...](#)

Paraffin Bäder

▶▶▶ [zur Kategorie...](#)

MoVeS Para Pro User Manual

1. About paraffin therapy

The MoVeS Para Pro Paraffin Bath is the perfect way to successfully apply thermal therapy for people with arthritis or other rheumatic diseases. The heat helps increase blood flow and relax the muscles, which can help relieve pain caused by arthritis, osteoarthritis and fibromyalgia. Paraffin wax can even soften hardened skin caused by scleroderma. It increases the skin's elasticity, allowing for increased movement and mobility, especially on the skin covering the hands. But the benefits of paraffin wax don't stop there. Because paraffin wax treatments are a form of heat therapy, they're often used for muscle, tendon and ligament ailments. As with conditions like arthritis, they increase blood flow, reduce joint stiffness and reduce pain. They're also used to treat bursitis, tendonitis, sprains and pulled muscles.

The MoVeS Para Pro Paraffin Bath can be used for aesthetic treatment as well, as in recent years has become very popular in Europe, America and South-East Asia. Its use is not only limited to the face, it can provide a whole body treatment. Paraffin treatment is the perfect method to make dry and rough skin become smooth and soft again, earning its nickname 'Crystal Mask' for the amazing results. The paraffin's vitamins and nutrients, which are released by melting, are beneficial for the skin.

When the melted paraffin comes into contact with the skin, it immediately solidifies, forming an insulating layer. This will help stimulate the circulation of blood, enhance nutrition absorption, soften wrinkles and moisten the skin, making it look radiant again. Paraffin treatment can easily be applied to hand, neck, chest, legs and feet.

2. To start

This paraffin bath is equipped with an electronic thermostat with a temperature range from 40°C to 60°C. Plug the unit into a grounded outlet of the correct voltage (220V, 50/60Hz). Be sure it's on a level surface and away from a sunny window, radiator or any other heat source.

Place the paraffin in the bath. To speed up the melting process, put the thermostat on "10" for one hour and then lower it to 56°C or less. Never insert any part of the body when the temperature is above 56°C!

This unit is designed to stay on 24 hours a day. The indicator light will go on and off intermittently during the normal operation. When not in use, always leave the lid on.

3. The treatment

Hand

1. Remove all jewelry. Thoroughly sanitize the skin.
2. Use tissue to completely dry the hand.
3. Apply a collagen emulsion on the hands and rub until absorbed.
4. Test the temperature before dipping. Gently dip the hand into the melted wax and pull it out. Repeat this action about 5 times, until a 5 to 10mm layer of wax has formed.
5. Put the plastic glove on the hand.
6. Lay a towel glove over the hand to retain the warmth.
7. Relax and leave the wax on for 10 minutes.
8. Peel off like a glove and discard.
9. Use a moisturizer immediately after removing the wax because the pores are open and able to absorb it more readily.

Face

1. Place a brush in the bath for five minutes to soften.
2. Clean the face, wipe it dry, then apply a collagen emulsion on the face and rub until absorbed.
3. Test the temperature, cover the eyes with cotton cloth.
4. Use the brush to apply the wax on face.
5. Lay a gauze on the face.
6. Apply the wax on the gauze for 3 to 4 times to achieve the best results.
7. Relax and leave the wax on for 15 minutes.
8. Peel off the gauze from the face and discard.
9. Use a moisturizer immediately after removing the wax because the pores are open and able to absorb it more readily.

Feet

1. Clean the foot, thoroughly sanitize the skin, wipe it dry, then apply a collagen emulsion on the foot and rub until absorbed.
2. Test the temperature before dipping. Gently dip the foot into the molten wax and pull it out. Repeat this action about 5 times, until a 5 to 10mm layer of wax has formed.
3. Put the plastic glove on the foot.
4. Lay a towel glove over the foot to retain the warmth.
5. Relax and leave the wax on for 10 minutes.
6. Peel off like a glove and discard.
7. Use a moisturizer immediately after removing the wax because the pores are open and able to absorb it more readily.

Note:

- For health consideration, use different paraffin baths for hands, feet and body treatment.
- Do not re-use wax.

4. Contraindications

- Do not use paraffin treatment in case of:
 - Open wounds
 - Skin irritation
 - Skin infection
 - Reduced thermal sensation.

5. Maintenance

Each time before use, check if the paraffin level is between the minimum and maximum indicator lines. When the level becomes too low, refill with paraffin to the correct level.

Do not turn on the power without paraffin in the bath.

Do not add other products to the paraffin (bath), as this may change the melting point and might affect the machine's command and control application.

6. Cleaning

- After 50 treatments, dirt and other undesirable impurities may collect at the bottom of the wax tank. When this happens, the wax should be removed.
- Cut the power off, wait until the paraffin becomes cool, dry and hard.
- To remove the solidified paraffin block from the unit, turn it on for a few minutes to 10 until the edges of the paraffin start melting and the sides loosen. Push one end of the block, making it tilt up. Lift the block out of the bath.
- Turn the power off and use a facial tissue to clean the bath. Use wax purifier to allow layers of impure wax to be removed thoroughly.
- Refill with new paraffin and start to use.

7. Warning

- Do not fill the paraffin bath with water or do not soak the paraffin bath in water. Doing so may cause damage beyond repair.
- In the event of a malfunctioning machine, please do not attempt to fix the machine yourself. Please contact your local reseller for check-up and repair.
- Please assure that the machine is on a dry and stable surface.
- Always use the paraffin bath in a well ventilated area.
- When not in use, always cover the tank with the lid.

8. Accessories detail

Purifier liquid, paraffin, cling film hands glove, cling film feet glove, towel, brush, gauze, collagen emulsion etc. mentioned in this manual must be purchased separately.

9. Specifications

Dimension	39cm x 25cm x 18cm (LxWxH)
Net weight	3.1 kg
Power supply	220V - 50Hz/60Hz
Power consumption	150W
Operation temperature	40°C – 60°C
Temperature accuracy	setting temperature $\pm 1.5^\circ\text{C}$
Paraffin Capacity	2 to 3 kg (min – max indicator level)
Operation time	24/7

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